



Storage up to 30°C



18 months


### Sample preparation for LM detection

#### Detection in solid/ liquid foodstuff and water

Add 9 times volume of half-Fraser broth to the sample and homogenize. Incubate at 30 +/- 1 °C for 25 +/- 1 hours to enrich the culture.

#### Viable count in swab test specimen

Add 1mL of wiping solution to 9 mL of half-Fraser broth. Incubate at 30 +/- 1°C for 25 +/- 1 hours to enrich the culture.



### CompactDry Protocol

**1** Open the cap **2** dispense 1 mL of sterile diluent (e.g. saline) on the middle of the CompactDry plate. **3** Pipette 0.1 mL of enriched culture in the middle of the plate. **4** To achieve single colonies, streak the inoculum gently from top to bottom using a loop. Spread it over the entire plate surface. **5** Replace the cap, label and turn the plate over. **6** Incubate for the required time and temperature. Following incubation, count the number of coloured colonies.

### Preparation for enumeration

#### Viable count in water or liquid foodstuff


1 mL of sample (diluted if necessary) deposited on the LM CompactDry plate.

#### Viable count in solid food

Add 9 part of buffered peptone water to 1 part of food sample and homogenize by stomacher®. 1 mL of sample (diluted if necessary) be dispensed in the middle on the LM CompactDry plate.

#### Viable count in swab test specimen

Use the swab to wipe the surface, put the swab back into the tube and screw tightly. Shake the tube to mix. 1 mL of sample (diluted if necessary) is dispensed in the middle on the LM CompactDry plate.



**Incubation temperature 37 +/- 1 °C**

**Incubation time 24 +/- 2 hours plus a further 24 +/- 2 hours**

#### Interpretation-Detection/Enumeration

Listeria monocytogenes form red colonies with or without blue halo. If colonies of L. monocytogenes are observed, perform confirmation tests by ISO11290-1:2017, ISO11290-2:2017 or other methods.

#### General information

Listeria ivanovii also forms red colonies with or without blue halo. Other Listeria spp. form blue/green colonies.

Bacteria other than Listeria spp. are inhibited by selective agents in the medium and do not form coloured colonies if they do grow. Rarely some Bacillus spp. may form relatively large, flat/orange colonies.

The growth area is 20 cm². Use the 1 cm² grid to ease colony counting.

High concentrations cause the entire growth area to become red coloured. In this case dilute the sample.

After use please follow local disposal regulations.

#### Storage and shelf life

Keep at room temperature (+ 1 to +30 °C). Shelf life 18 months after manufacturing.